



CENTRAL MEDICAL SERVICES SOCIETY

(An Autonomous Body under Ministry of Health & Family Welfare, Govt. of India)

2nd Floor, Vishwa Yuvak Kendra, Teen Murti Marg,
Chanakypuri, New Delhi-110021.



Date: - 09/12/2025

Minutes of Pre-bid meeting held on 27/11/2025 at 11:00 AM regarding Tender No.: CMSS/PROC/2025-26/FP/025, GeM Bid No: GEM/2025/B/6917249, dated 21/11/2025 for Procurement of Social Marketing Contraceptives for Family Planning

Pre-bid queries raised by the prospective bidders and Clarification/ Amendments by CMSS:

Sr. No.	As per tender	Tender clause no. & Page No.	Bidder's Query	Bidder's Name	Clarifications/ Amendments
1.	Technical Specification of Schedule I & II	Section VI & Pg No. 59	Foil Specification: As per Tender document, the specification is mentioned as Paper-based foil, whereas our last supplies were made with PET-based foil Refer attached Annexure A # 1 (Page No.59) - Tender specification Kindly specify the foil specification for Schedule I & II.	M/s TTK Healthcare Limited	Clarified as: Kindly adhere to the technical specification mentioned in the tender.
2.			Foil Width: As per Tender document, the specification mentioned as 70+2mm length x maximum 30mm width, whereas our last supplies were made with 55mm length x 55mm width Refer attached Annexure A # 2 (Page No.59) - Tender specification Kindly specify the foil width for Schedule I & II.		Clarified as: Kindly adhere to the technical specification mentioned in the artwork. Square pack of strips dimension: 55 x 55 mm (L x B) Quantity of Condom: 1 pc. condom/ square strip
3.			Wallet Dimension: As per Tender document, the specification mentioned as 75 x 50 x 13mm, whereas our last supplies were made with 65 x 25 x 65mm Annexure B: Attached artwork of last supplies - Schedule I. Refer attached Annexure A # 3 (Page No.59) - Tender Specification		Clarified as: Kindly adhere to the technical specification mentioned in the artwork. Dimension - 58 x 60 x 23 mm (LxBxH) Quantity of Condom: 5 square strips/wallet
4.			Carton Dimension: As per Tender document, the specification mentioned as		Clarified as:

Tel. No. 011-21410905/6, Fax No. 011-21410849 Website: www.cmss.gov.in



		145 x 120 x 80mm (20 wallets), whereas our last supplies were made with 205 x 135 x 110mm (24 wallets x 5 pcs = 120 pcs) Refer attached Annexure C # 1 (Page No.60) - Tender Specification Refer attached Annexure D - our last supplies Carton specification	Kindly adhere to the technical specification mentioned in the artwork. Dimension - 182 x 195 x 127mm Quantity of Condom: 5 x 48 = 240 pcs.
5.		Outer Box Dimension: As per Tender document, the specification mentioned as 490 x 300 x 410mm, whereas our last supplies were made with 570 x 425 x 425mm (30 outer x 24 wallet x 5 pcs = 3600 pcs in each outer box) Refer attached Annexure C# 2 (Page No.60) - Tender Specification Refer attached Annexure E- our last supplies Outer Box specification	Clarified as: Kindly adhere to the technical specification mentioned in the artwork. Dimension - 555 x 408 x 525mm Quantity of Condom: 5 x 24 x 45 = 5400 pcs.
6.		Outer Box Quantity: As per Tender document, the specification mentioned as 40 cartons of 30 wallets x 5 nos. of Deluxe Nirodh (i.e. 4000 pcs) whereas the actual calculation of 40 cartons of 30 wallets x 5 nos. comes to 6000 pcs. Refer attached Annexure C # 3 (Page No.60) - Tender Specification	Quantity of Condom: 5 x 24 x 48 = 5760 pcs.

Note: - Above changes will be part of the tender document and apart from above, all other terms & conditions of tender document will remain unchanged.


 (Manish Kumar Bhaskar)
 Joint Director (Procurement)





CENTRAL MEDICAL SERVICES SOCIETY

(An Autonomous Body under Ministry of Health & Family Welfare, Govt. of India)

2nd Floor, Vishwa Yuvak Kendra, Teen Murti Marg,
Chanakyapuri, New Delhi-110021.



Date: - 09/12/2025

CORRIGENDUM No.-1

Tender for Procurement of Social Marketing Contraceptives for Family Planning

Tender No.: CMSS/PROC/2025-26/FP/025

GeM Bid No: GEM/2025/B/6917249, dated 21/11/2025

1. The Tender timeline for the above referred tender has been amended as follows:

Particulars	Existing	Amended as
Last date and time for bid submission	15/12/2025 at 03:00 PM	16/12/2025 at 03:00 PM
Last date and time for submission of original documents	15/12/2025 at 03:00 PM	16/12/2025 at 03:00 PM
Date and time for tender opening (technical bid)	15/12/2025 at 03:30 PM	16/12/2025 at 03:30 PM

Note: - Apart from above, all other terms & conditions of tender document will remain unchanged.


(Manish Kumar Bhaskar)
Joint Director (Procurement)



2. The artwork of Schedule I & II is attached herewith.

Condoms Deluxe Nirodh (HLL)

DELUXE NIRODH PACKING SPECIFICATIONS

FOIL:

STRIP : 4 Ply Laminated Foil

Material : 40/42 gsm GIP(glassine) paper/0.009mm alluminium foil/ heat scal coating/
Acid Copolymer of polythylene

B Type : Strip size for non-squeeze pack, minimum 55mm length x 55mm width.

Printing : In 2 colours, black and red on glassine as per existing design.



WALLETS - box type

- Quantity : (5 condoms) directly placed in wallet.
- Material : 250 gsm duplex board printed in four colours (as per existing design).
- B Type : Size of the wallet minimum dimensions 58 x 23 x 60mm
- Printing : Printing on wallets in four colours on off set press. The wallets finishing of double coat scuff proof, varnish cut and creased. The off set printing will be done based on the artwork supplied.





INSTRUCTIONS - Leaflet type

Quantity : 1 Leaflet directly placed in wallet.

Printing : In 2 colours, Black & White (as per existing design).

5.5cm x 5.5cm



CARTON – shelf box type

Quantity : 48 wallets (240 condoms) per shelf box.

Material : 3 ply flute 100 x 100 x 250 duplex board top open.

B Type : Size of the carton minimum dimensions 182 x 195 x 127mm.

Printing : In 4 colours (as per existing design).




OUTER BOX - box type

Quantity : 24 shelf boxes (5760 condoms) per outer box.

Material : Narrow flute corrugated card board boxes of 7 ply each ply of 150 gsm virgin kraft paper of which outer ply will be alkali resistant.

B Type : Size of the carton minimum dimensions 555 x 408 x 525mm.

Printing : In 2 colours, Red & Black (as per existing design).



Condoms Josh (HLL)



Integrity Packaging LLP

911, 9th Floor, Gold Crest Building, Opp. Tanishq Jewellery,
L.T. Road, Borivali (W), Mumbai 400 082, Maharashtra, India.
T: +91-22-43471770 - 1771
Email: service@integrityprint.in | Web: www.integrityprint.in

Client : **HLL Lifecare Limited (Chennai)**


Date: **18.05.24**


Item Name : **IPLLP-01431 Josh Wallet 3 Pcs**

Printing Colour Code




Reference Colour Code





Job Name :- Josh 3s Tray
Job Size :- L120 X W 40X H 235
Date :- 09-04-2024



Josh 3s
235X120X40 mm (LxBxH)

Condoms Masti (PHSI)



Back Side Printing

IMPORTANT

1. This product is for single use only. Use a new condom every time you have sex & follow the instructions carefully.
2. Make sure there is an adequate lubrication. If you add lubricant, use a water based lubricant. DO NOT USE OIL BASED LUBRICANTS such as those made with petroleum jelly (e.g., Vaseline), mineral oil, vegetable oil or cold cream as they may damage the condom.
3. "Condom greatly reduces the risk of unintended pregnancy, HIV/AIDS and sexually transmitted diseases if used correctly during sex."
4. Do not use two condoms together during sex because condoms will rub against each other and may get tearing or break.
5. Seek medical assistance as soon as possible, at least within 72 hrs. if the condom leaks/breaks during use.
6. If you are using any medical treatment in the genital area (creams, lotions, ointments etc.), consult your doctor to ensure that the treatment won't affect the condom adversely.
7. Do not use condom after expiry date, indicated on the pack.
8. Any use of latex condom for other than vaginal intercourse may increase the potential of damage to the condom.
9. Do not reuse the condom, if reused there might be an increased risk of failure or infection.
10. The condom is made of natural rubber latex which might cause allergic reactions including anaphylactic shock.
11. This product complies with all requirements of ISO 4074.
12. Please stop and check if the user feels the condom slipping off or tightening excessively on the penis because the right size is crucial.

INSTRUCTIONS FOR USE

- 1) Use a new condom for each act of sex. Check the condom if pack is found damaged.
- 2) Tear to open the foil carefully. Do not use fingernails, teeth or anything that can damage the condom.
- 3) Before any sexual act, place the condom on the head of erect penis ensuring the tip of the condom is inside.
- 4) Squeeze the tip of condom to remove air. Excess air could cause the condom to burst. Latexes, use-squeezable condoms, semen, vaginal secretions and blood can all transmit infectious organisms.
- 5) Now slowly unroll it from the tip of the head to the body of the penis. If the condom does not unroll easily, it may be damaged or too old. Do time it out and start with a new condom.
- 6) Immediately after ejaculation, hold the rim of condom and withdraw from another point.
- 7) Remove the used condom by wrapping it in a tissue and throw it into the dustbin. Wash your hands properly with soap and water.
- 8) Do not throw the used condom into the Toilet to avoid any blockage.



DOI: 10.1002/for

139 mm


184 mm

195 mm




Condoms Rakshak (HLL)

Job Name :- Rakshak Wallet 2s 08
Job Size :- L 58 X W 10 X H 60
Date :- 10-04-2024



AW No.HLLPFT/WGR2/2101/Rev.02/DEC.2022

Job Name :- Rakshak 2S Try
Job Size :- L 120 X W 40 X H 260
Date :- 09-04-2024



Rakshak 2s
260X120X40 mm (LxBxH)

Condoms Sawan 3s & 10s (PSS)

ANNEXURE-1

006/02d/0624/3PCS




उपपर सचिव/Under Secretary
स्वास्थ्य एवं परिवार कल्याण मंत्रालय
Ministry of Health & F.W.
भारत सरकार/Govt. of India
नई दिल्ली/New Delhi

ANNEXURE-2

006/02d/0624/10PCS

Sawan New 03 Bd
10 Pcs



10/7/24
 (NOUN TOPNO)
 अवर सचिव / Under Secretary
 स्वास्थ्य एवं परिवार कल्याण मंत्रालय
 Ministry of Health & F.W.
 भारत सरकार / Govt. of India
 नई दिल्ली / New Delhi

Condoms Thrill (PHSI)





SMOOTH
GELDY ZINCHET BENTHRA WOODS
THRIL
®

60 x 3 Pcs



THRIL
SMOOTH PREMIUM LUBRICATED
CONDOMS

60 x 3 Pcs

THRIL
SMOOTH PREMIUM LUBRICATED
CONDOMS

DIRECTIONS FOR USE

1. Carefully open the Pack, Unwrap Pen or out of the tip
2. Slowly unroll the Condom on the base Pen. This is essential before any Sexual Entry. 3. After withdrawal, slide the condom off without springing the condom without allowing the Penis to touch vagina
4. Use condoms to prevent pregnancy, STD and HIV/AIDS infection, gonorrhoea, chlamydia, syphilis, hepatitis B, HIV, AIDS and other sexually transmitted diseases.

STORAGE: IN A COOL DRY PLACE AWAY FROM HEAT,
DIRECT SUNLIGHT OIL AND MECHANICAL DAMAGES

Subsidiary by Govt. Of India under GMP
Not for sale outside India
Pack : Type A - Square Condom Length : 170 mm (Min.)
Nominal width : 48 mm ± 2 mm

Pink Lubricated Electrostatically Tested Condoms
For Single Use Only



THRIL
SMOOTH PREMIUM LUBRICATED
CONDOMS

60 x 3 Pcs

थ्रिल
स्मूथ प्रीमियम लुब्रिकेटेड
कांडम

60 x 3 Pcs

Manufactured by:
Indus Medicare Private Limited
Ramapally (V), Medak District,
Telangana State, India - 502336
Mfg. Lic.No. BPG/MO/2023/04293

Marketed by:
Population Health Services (India)
2A, SOA Compound, Jawahar Mera
Dewas Naka Indore (Madhya Pradesh)

Consumer Call : Admin Manager
Post Box No. 1523, Hyderabad - 500 062
Phone : 040-25166044
e-mail : feedback@pshindia.org.in

Batch No.
Mfg. Date,
Expiry Date,
Lot No.



MRP ₹ 10/- for 3 Pcs
(Inclusive of all taxes)

अधिकतम सुरक्षा के लिए कांडम को एक ही बार में ही प्रयोग करें।
USP (Inclusive of all taxes) ₹ 3.34 Per Condom
थ्रिल कांडम का उपयोग एक ही बार में ही करें।

Condoms Ustad (HLL)

 Integrity Packaging LLP 311, 3rd Floor, Gold Crest Building, Opp. Taramani Jewellery, L.T. Road, Survali, TN, Mumbai 400 092, Maharashtra, India. T: +91 22 45411792 - 1771 Email: dmw@iprint.in Web: www.integrityprint.in	Client : HLL Lifecare Limited (Chennai)		Date: 05.04.24																								
	Item Name : Ustad Wallet 3 Pcs																										
	<table border="1"> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td>Cyan</td> <td>Magenta</td> <td>Yellow</td> <td>Black</td> </tr> </table>										Cyan	Magenta	Yellow	Black	Printing Colour Code <table border="1"> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </table>												
		Cyan	Magenta	Yellow	Black																						
		Reference Colour Code <table border="1"> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </table>																									





Integrity Packaging LLP

911, 9th Floor, Gold Crest Building, Opp. Tanishq Jewellery,
L.T. Road, Borivli (W), Mumbai 400 092, Maharashtra, India.
T: +91 22 43471770 - 1771
Email: artwork@integrityprint.in | Web: www.integrityprint.in

Client : HLL Lifecare Limited (Chennai)


Date: 05.04.2024

Item Name : Ustad Wallet 6 Pcs

Printing Colour Code



Reference Colour Code





[illegible]

- 58mm

15mm

65mm



58mm

14mm

Zaroor
CONDOMS

Zaroor condoms are pink coloured, lubricated, top-ended, 100% electrically tested, natural rubber latex condoms with a vanilla masking scent. 53 ± 2 mm nominal width and 190 mm (min.) length. Each condom is meant for single use only.

Zaroor condoms are manufactured to meet the highest international Quality Standards to provide reliable protection.

Zaroor is a brand of DKT, which is based in USA and has affiliate offices in Latin America, Middle East and Asia.

Marketed by: DKT India (An Affiliate of DKT International, USA)
Hem Dill 67-A, Linking Road, Santacruz (West), Mumbai 400 054.
Customer Care, Manager - Logistics, Tel. No: 8455954776
Email: customercare@dktdindia.org

Manufactured by: HSA Lifecare Limited
(A Government of India Enterprise) Perambur,
Thiruvananthapuram 695 905, Kerala.

Mfg. Lic. No.: MF6/MS/2822/988644

Subsidized by Government of India under CSMP, Central Government Safety Net For Export Goods India.

MRP (incl. of all taxes)
shown above per
pair (ex. sales)

₹ 20/-

Net Quantity : 6 Nos

Batch No. (see vial)

Mfg. Date (see vial)


Exp. Date (see vial)



LATEX

dkd
INDIA

65
PACK



Instructions for use:
• Do not tear or puncture the condom before use.
• If not erect, roll the condom back and unroll the tip of the condom and unroll.
• After sex, hold the condom in place and withdraw the penis while it is still erect.
• Tie a knot and throw the condom in the trash.

113mm

Zaroor condoms are sufficiently lubricated. Use only water based lubricants if required. Oil based lubricants like cooking oil, Vaseline may cause damage to the condom.



